



- DE** Anwendungen
- EN** Application
- FR** Application
- ES** Aplicación



DE	Polierer	Diamantpolierer
EN	Polishers	Diamond Polishers
FR	Polissoirs	Polissoirs Diamantés
ES	Pulidores	Pulidor Diamantado

9825	9826	9827	9828	9831	9832	9833	9834			
9812G	9813G	9814G	9812M	9813M	9814M	9812F	9813F	9814F		
9501G	9502G	9504G	9501M	9503M	9504M	9507M	9501F	9504F	9144F	
9140EF	9141EF	9142EF	9143EF	9145EF						
9109M	9111M	9112M	9113M	9114M	9115M	9119M				
9133G	9134G	9133M	9134M	9133F	9134F					
9121M	9122M	9123M	9124M	9121F	9122F	9123F	9124F			
ZIR9861M	ZIR9862M	ZIR9861F	ZIR9862F	ZIR9863M	ZIR9863F					
LS9871M	LS9872M	LS9871F	LS9872F	LS9873M	LS9873F					
9150										
9837										





	✓	✓	✗
DE	empfohlen	geeignet	nicht geeignet
EN	recommended	suitable	not suitable
FR	recommandée	adaptée	ne convient pas
ES	recomendado	posible	no es posible



Keramik	Zirkonoxid	Lithium-Disilikat	Komposit	Gold / Amalgam	Chrom Legierungen / Edelstahl	Schmelz	Kunststoff	Nano Keramik
Ceramic	Zirconium oxide	Lithium disilicate	Composite	Gold / Amalgam	Chromium precious / stainless steel	Enamel	Plastic	Nano Ceramic
Céramique	Dioxyde de zirconio	Di silicate de lithium	Composite	Or / Amalgame	Chrome précieux / acier inoxy-dable	Émail des dents	Matière plastique	Nano céramique
Cerámica	Dióxido de zicóno	Di silicado de litio	Composite	Oro / Amalgama	Cromo preciosos / acero inoxidable	Esmalte	Plástico	

✗	✗	✗	✓	✗	✗	✗	✗	✗
✓	✓	✓	✗	✗	✓	✗	✗	✗
✗	✗	✗	✗	✓	✓	✗	✗	✗
✗	✗	✗	✓	✓	✗	✗	✓	✗
✓	✗	✗	✓	✓	✓	✓	✓	✗
✓	✗	✗	✗	✗	✗	✗	✗	✗
✗	✗	✗	✓	✗	✗	✗	✗	✗
✗	✓	✗	✗	✗	✗	✗	✗	✗
✓	✗	✓	✗	✗	✗	✗	✗	✗
✗	✗	✗	✗	✓	✓	✗	✗	✓
✗	✗	✗	✓	✗	✗	✗	✗	✓

