

jota podiatry
Safe-Edge Instruments

**New tools for accurate and injury-free
foot care & foot treatment!**

- › „Safe Edge“ Design prevents from injuries -
Perfect for risk patients (diabetics)
- › Stainless Steel, Natural Diamonds and Swiss
Precision - for successful treatment & superior
results!



jota podiatry Safe-Edge Instruments

“Safe Edge” design provides utmost safety during therapies and prevents skin and nail bed from injuries.

Top Grinder P837T.104.052 806104179523052

For a selective treatment of diseased, thickened or deformed nails:

- On top diamond coating only
- Special rounded-edge pass from top to side
- No contact of grinding surface with adjacent skin
- Reduced risk for injuries
- Ideal for diabetics patients



Side Grinder P837S.104.055 806104147523055

Ideal for a safe therapy at nail borders and close to nail beds:

- Smooth, No-grinding Top
- Special rounded-edge pass from top to side
- Safe treatment at border regions to skin and nail bed
- No undesired skin contact
- Reduced risk for injuries
- Ideal for diabetics patients



Smooth, round edge,
not damaging the skin

Angular edge,
may irritate, damage the skin



Round-Edge Grinder P837R.104.055 806104157523055

The perfect addition for skin removal, grinding and smoothing:

- Round edge on pass from top to side
- Safe treatment
- Perfect for scarf skin and nail surfaces



Recommended Rotation speed: 5.000 - 10.000 min⁻¹

We recommend a decrease in rotation speed during dry use and or risk-patients in order to prevent from temperature increase.